

DNA Ladder 25 - 700 bp




Cat. N°: DG6001 (50µg) / DG6001-XL (250µg)

Description:

Low Range Ladder is a ready-to-use DNA Marker suitable for sizing linear double-stranded DNA fragments from 25-700bp, and composed of 10 linear chromatography-purified individual DNA fragments. All bands (except 100bp and 300bp) are supplied at approximately 40ng/5µl. The 100bp and 300bp bands have increased intensity relative to the other bands on an ethidium bromide-stained agarose gel.

Concentration: 0.1 µg/µL

Storage condition: Stable at room temperature. For long term storage, store at -20°C.



PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development and is not suitable for administration to humans or animals.