

## DNA Ladder 50 - 1,000 bp




Cat. Nº: DG6002 (50µg) / DG6002-XL (250µg)

### Description:

The 50-100 bp DNA ladder is suitable for sizing double stranded DNA from 50 bp to 1kb. The DNA ladder is ready to use and can be visualized on agarose gel. Consist of 14 fragments: One step ladder ranging from 50-500 base pairs and another ranging from 100-500 base pairs. The resulting steps are picture in the figure. **Concentration:** 0.1 µg/µL

**Storage condition:** Stable at room temperature. For long term storage, store at -20°C.



### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development and is not suitable for administration to humans or animals.