


## Lambda DNA/HindIII Marker



Cat. Nº: DG6007 (50µg) / DG6007-XL (250µg)

### Description:

The Lambda DNA/HindIII marker is ready to use and can be visualized on agarose gel. It is composed by 8 fragments ranging from 125 to 23.13 kilobases in size: 125, 564, 2027, 2322, 4361, 6557, 9416, 23130 bp. The resulting steps are shown in the figure

**Concentration:** 0.1 µg/µL

**Storage condition:** Stable at room temperature. For long term storage, store at -20°C.



### PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development and is not suitable for administration to humans or animals.